

## Immunoprecipitation of fusion proteins from cell extracts using NanoTag's Selector Resins

NanoTag's Selector Resins are based on high-affinity single-domain antibody (sdAb) fragments covalently immobilized on 4 % cross-linked agarose beads. The innovative, oriented and selective attachment via flexible linkers guarantees optimal accessibility of the sdAbs and in addition largely eliminates batch-to-batch variations. Due to the single-chain nature of sdAbs and their stable and covalent attachment, no leakage of light and heavy chains is observed during elution with SDS sample buffer. Selector Resins thus feature high affinity and superior capacity for fusion proteins while showing negligible non-specific background.

Selector Resins are compatible not only with physiological buffers but also with high stringency buffers (see product-specific compatibility charts) and reducing agents. Selector Resins thus provide great freedom to adjust the binding and washing conditions to the experimental needs.

**Only for research applications, not for diagnostic or therapeutic use!**

### Products

<b>MBP Selector</b>	N0110	2000 µL slurry (1000 µL resin)
<b>GST Selector</b>	N0210	2000 µL slurry (1000 µL resin)
<b>GFP Selector</b>	N0310	2000 µL slurry (1000 µL resin)
<b>RFP Selector</b>	N0410	2000 µL slurry (1000 µL resin)
<b>TagFP Selector</b>	N0510	2000 µL slurry (1000 µL resin)

### Materials (not included):

- Lysis buffer<sup>(1)</sup>
- Washing buffer<sup>(1)</sup>
- Tris-buffered saline (TBS) pH 7.4
- 2x SDS sample buffer

### Remarks

<sup>(1)</sup> Selector resins are compatible with most common Lysis/Washing buffers (e.g. RIPA). For custom buffers please refer to the product-specific compatibility chart.

**For further information concerning this protocol please contact us at  
[info@nano-tag.com](mailto:info@nano-tag.com)**

---

## Protocol using Mini Spin Columns

---

**We recommend using Mini Spin Columns (Cat. No. A1001) for immunoprecipitation experiments. An alternative batch protocol is found on the following page.**

1. Prepare native cell lysates (0.2 to 1.5 ml volume) according to established protocols. For mammalian cells, we recommend using  $10^6$ - $10^8$  cells per experiment.
2. Clarify lysate by centrifugation C for 10 min at  $> 14000 \times g$  and  $4 \text{ }^\circ\text{C}$ .  
Take sample for further analysis (Input fraction)
3. Equilibrate GFP Selector Resin
  - a. Resuspend GFP Selector Resin.
  - b. Transfer  $20 \text{ }\mu\text{l}$  slurry ( $10 \text{ }\mu\text{l}$  packed beads) into a clean 1.5 ml reaction tube.
  - c. Add 1 ml lysis buffer.
  - d. Centrifuge for 1 min at  $1000 \times g$  and carefully remove supernatant.
  - e. Repeat steps c-d once.
4. Add clarified lysate from step 2 to equilibrated GFP Selector Resin obtained in step 3.
5. Incubate 1 h at  $4 \text{ }^\circ\text{C}$  with head-over-tail rotation.
6. Sediment beads by centrifugation for 1 min at  $1000 \times g$  and  $4 \text{ }^\circ\text{C}$ .  
Take sample from supernatant for further analysis (Non-bound fraction)
7. Washing
  - a. Carefully remove supernatant.
  - b. Resuspend beads in 1 ml Lysis buffer
  - c. Centrifuge for 1 min at  $1000 \times g$ .
  - d. Remove Supernatant
8. Transfer
  - a. Remove bottom plug from Mini Spin Column. Place column in 2 ml reaction tube.
  - b. Resuspend beads in  $200 \text{ }\mu\text{l}$  Lysis buffer, transfer suspension to Mini Spin column
  - c. Wash out beads sticking to tube with  $200 \text{ }\mu\text{l}$  Lysis buffer and transfer to column.
  - d. Centrifuge column for 1 min at  $1000 \times g$ , discard flow-through.
9. Wash twice with  $400 \text{ }\mu\text{l}$  Washing buffer, centrifuge for 1 min at  $1000 \times g$ .
10. Wash once with  $400 \text{ }\mu\text{l}$  TBS, centrifuge for 1 min at  $3000 \times g$ .
11. Attach bottom plug and place Mini Spin Column in a clean 1.5 ml reaction tube
12. Resuspend GFP Selector resin in  $50 \text{ }\mu\text{l}$  2x SDS sample buffer
13. Heat Mini Spin Column to  $95 \text{ }^\circ\text{C}$  for 2 min.
14. Remove bottom plug and centrifuge for 1 min at  $3000 \times g$ .  
Boil collected eluate for 5 min at  $95 \text{ }^\circ\text{C}$  and analyze by SDS-PAGE.

---

## Batch protocol

---

**Due to the more effective washing steps, we recommend using Mini Spin Columns (Cat. No. A1001) for immunoprecipitation experiments (see protocol on page 2).**

1. Prepare native cell lysates (0.2 to 1.5 ml volume) according to established protocols. For mammalian cells, we recommend using  $10^6$ - $10^8$  cells per experiment.
2. Clarify lysate by centrifugation for 10 min at  $> 14000 \times g$  and  $4 \text{ }^\circ\text{C}$ .  
Take sample for further analysis (Input fraction)
3. Equilibrate GFP Selector Resin
  - a. Resuspend GFP Selector Resin.
  - b. Transfer  $20 \text{ }\mu\text{l}$  slurry ( $10 \text{ }\mu\text{l}$  packed beads) into a clean 1.5 ml reaction tube.
  - c. Add 1 ml lysis buffer.
  - d. Centrifuge for 1 min at  $1000 \times g$  and carefully remove supernatant.
  - e. Repeat steps c-d once.
4. Add clarified lysate from step 2 to equilibrated GFP Selector Resin obtained in step 3.
5. Incubate 1 h at  $4 \text{ }^\circ\text{C}$  with head-over-tail rotation.
6. Sediment beads by centrifugation for 1 min at  $1000 \times g$  and  $4 \text{ }^\circ\text{C}$ .  
Take sample from supernatant for further analysis (Non-bound fraction)
7. Washing
  - a. Carefully remove supernatant.
  - b. Resuspend beads in 1 ml Lysis buffer
  - c. Centrifuge for 1 min at  $1000 \times g$ .
  - d. Wash beads 2-3 times with Washing buffer.
  - e. Wash beads once with TBS
8. Transfer beads in clean 1.5 ml reaction tube.
9. Centrifuge for 1 min at  $3000 \times g$ .
10. Carefully and completely remove supernatant
11. Resuspend GFP Selector resin in  $50 \text{ }\mu\text{l}$  2x SDS sample buffer
12. Heat for 5 min to  $95 \text{ }^\circ\text{C}$ .
13. Centrifuge for 1 min at  $3000 \times g$ .  
Collect supernatant and analyze by SDS-PAGE.